Whatman™ Polycap™ TF Disposable Filter Capsule

Product Information sheet

Warning

For research use only.

**Not recommended or intended
for diagnosis of disease in
humans or animals.**

**Do not use internally or
externally in humans or
animals.**

- Solvent filtration.
- Solvent resistant capsule.
- Hydrophobic PTFE membrane.
- Autoclavable.
- Lot & Catalog No. stamped into each housing.

These products have been designed and manufactured with PTFE membranes and polypropylene housings. No glue, adhesive, metal, epoxy or other extraneous materials are used in construction. All seals are fused.

Disposable filtration devices provide greater labor efficiency while ensuring superior filtration when compared to hand assembled reusable filter housings.

This document provides general information on the products listed. The specifications in the Technical Data section are intended to provide the basis for establishing functional use, as well as setting the quality assurance test performance levels.

Chemical resistance summary*

Classes of Substances 20°C (68°F)	Polypropylene/PTFE Guide* for use
Acetic Acid, 5%	R
Acetic Acid, Glacial	R
Acetone	R
Acetonitrile	R
Ammonia, 6N	R
Amyl Acetate	R
Amyl Alcohol	R
Benzene*	R
Benzyl Alcohol*	R
Boric Acid	R
Butyl Alcohol	R
Butyl Chloride*	R
Carbon Tetrachloride*	R
Chloroform*	R
Chlorobenzene*	R
Citric Acid	R
Cresol*	R
Cyclohexanone	R
Cyclohexane*	R
Diethyl Acetamide	R

Dimethyl Formamide	R
Dioxane	R
DMSO	R
Ethanol	R
Ethers*	R
Ethyl Acetate	R
Ethylene Glycol	R
Formaldehyde*	R
Freon TF*	R
Formic Acid	R
Hexane	R
Hydrochloric Acid, Conc*	R
Hydrofluoric Acid*	R
Isobutyl Alcohol	R
Isopropyl Alcohol	R
Methanol	R
Methyl Ethyl Ketone	R
Methylene Chloride*	R
Nitric Acid, Conc*	R
Nitric Acid, 6N*	R
Nitrobenzene*	R
Pentane*	R
Perchloro Ethylene*	R
Phenol 0.5%	R
Pyridine	R
Sodium Hydroxide, 6N	R
Sulfuric Acid, Conc*	R
Tetrahydrofuran*	R
Toluene*	R
Trichloroethane*	R
Trichloroethylene*	R
Water	R
Xylene*	R

R = Resistant; * = Short Term Resistance of Housing

*The above data is to be used as a guide only. Testing prior to application is recommended.

Typical Applications for Polycap TF Disposable Filter Capsules

Venting: Vent Filling, Mixing, Holding, Autoclaves and Fermentation Vessels.

In line Filtration: Sterilize and Remove Particles from Air/Gas/Liquid Streams. Remove Liquid Aerosols from Air/Gas Streams. Degas, Clarify and Sterilize Liquid Streams.

Other Uses: Clean Nitrogen used to dry samples. Ensure Clean Inert Gases in vial filling lines. Instrument Gas/Air preparation. Point of use "Ultra Clean" Gases, Air, solvents for Air Guns and Cleaning Wands.

Isolation (gases pass, liquids/aerosols stopped): Isolation of Gases or Vents for Incubators, Autoclaves, Lyophilizers, EtO Sterilizers. Isolation of Exhaust from Fermentors

Electronics: Photoresists, Solvents, Gases for Research or Industry.

Pharmaceutical: Vents and In Line for Production Systems.

Biotech: Vents and Exhausts for Biological Environments, in-line Sterile Gases.

Operating Instructions: Venting

Safety: Considering the special factors of your application, consult the table of Technical Data to determine the correctness of use. Key safety concern is to not exceed the pressure, temperature, or chemical compatibility recommendations.

PTFE Membrane Considerations: PTFE membrane is hydrophobic and will not allow water (aqueous solutions) to pass without high pressure. This pressure is called the Water Breakthrough Value (WBT) Value, and changes with the pore size of the membrane. Aqueous solutions may be filtered if the membrane is initially "wetted" with alcohol or other appropriate solvents. PTFE membrane will stop aqueous aerosols in gas streams.

Efficiency: To maximise filtration throughput, use the largest pore size

filter that will provide the required cleanliness. (For sterilization of liquids 0.2 µm is required.) If the material to be filtered is heavily contaminated, consider prefilters. To extend filter life use low flow or pressure.

Autoclaving: Autoclave; at 121°C (132°C Max.) for 20 minutes. May tolerate 132°C for one hour. Multiple autoclave cycles may be possible, however, reuse is the responsibility of the operator, who should protect the device from cross contamination and detect loss of integrity by appropriate testing. (EtO sterilization should be possible but no studies have been conducted. Radiation is destructive to PTFE.)

Filtration Installations:


Vents: We recommend the use of PolyVENT™ for sterile vent applications. Attach one connector to vessel, the other connector is open to the atmosphere. If exhaust gas is saturated with moisture, install vent filter in a vertical position to allow collected moisture to drain back into vessel rather than fill (and block) vent filter housing. Sterilizing vents are typically 0.45 µm, or for longer duration and maximum security, use a 0.2 µm membrane.

Air/Gas "In-Line": Effective service life varies. Replacement is required on a timely (routine) basis. Monitor pressure buildup or flow restriction to determine frequency. Pore size selection is based on application.

Solutions: Small volumes may be filtered using a syringe; larger volumes may be filtered by connecting to the outlet of a pressure vessel or "in line" using gravity, vacuum, or peristaltic pump for the pressure required to filter. **Air Locks** seriously hamper flow rates; to eliminate during initiation of flow, use low pressure (<5 psi/0.3 bar), bleed off air with vent, and then increase pressure and flow. The air in the system and filter device must be flushed prior to the membrane being wetted. (Use the manual vent to exhaust entrained air.)

Integrity Testing:

Bubble Point (BP) Test: Flush the PTFE membrane filter device with an appropriate solution (methanol). After the membrane is completely wetted, apply air under controlled pressure until air breaks through the membrane and bubbles from the exit. The pressure at which air passes through the "wetted" membrane is the "bubble point". Refer to the table for pore size and methanol BP values. The test is best conducted with outlet connector pointed upward.

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Design Features

- **0.1 µm** filter devices for "ultra clean" applications
- **Barbed Hose Connections** for multiple tubing sizes
- **Integrity testable: Bubble Point**
- **Vent on inlet**

Maximum Pressure: 4.1 bar (60 psi)



Housing: Polypropylene

Filter Media: PTFE

- **Compact:** Incredibly efficient design provides high effective filtration area in small size.
- **Package Clean "Double Wrapped"**
- **Manufactured in Clean Room GMP Controlled production facilities.**

Autoclaving: Autoclave at 121°C (132°C Max.) for 20 minutes. Multiple autoclave cycles may be possible, however reuse is the responsibility of the operator, who should protect the device from cross contamination and detect loss of integrity by appropriate testing.

Typical Flow Rates, Polycap TF

Connectors

● 1/4-3/8, 1/4-3/8

▲ 3/8-1/2, 3/8-1/2

✕ 1/4 MNPT, 1/4-3/8

■ 1/4 MNPT, 3/8-1/2

Catalog Number	Product	FILTER			Connections (1)		Methanol Bubble Point (2)		Qty/pkg
		Rated μm	Media	Area	Inlet	Outlet	Bar	psi	
6700-3602	Polycap 36 TF	0.2	PTFE	500 cm ² (78 in ²)	6-10 mm (1/4-3/8") SB	6-10 mm (1/4-3/8") SB	0.9	13	1
6700-3610	Polycap 36 TF	1.0	PTFE	500 cm ² (78 in ²)	6-10 mm (1/4-3/8") SB	6-10 mm (1/4-3/8") SB	0.2	3	1
6710-3602	Polycap 36 TF	0.2	PTFE	500 cm ² (78 in ²)	10-12 mm (1/4-3/8") SB	10-12 mm (3/8-1/2") SB	0.9	13	1
6711-3601	Polycap 36 TF	0.1	PTFE	500 cm ² (78 in ²)	1/4" MNPT	6-10 mm (1/4-3/8") SB	1.7	25	1
6711-3602	Polycap 36 TF	0.2	PTFE	500 cm ² (78 in ²)	1/4" MNPT	6-10 mm (1/4-3/8") SB	0.9	13	1
6711-3604	Polycap 36 TF	0.45	PTFE	500 cm ² (78 in ²)	1/4" MNPT	6-10 mm (1/4-3/8") SB	0.5	7	1
6700-7502	Polycap 75 TF	0.2	PTFE	1000 cm ² (155 in ²)	6-10 mm (1/4-3/8") SB	6-10 mm (1/4-3/8") SB	0.9	13	1
6700-7504	Polycap 75 TF	0.45	PTFE	1000 cm ² (155 in ²)	6-10 mm (1/4-3/8") SB	6-10 mm (1/4-3/8") SB	0.5	7	1
6701-7510	Polycap 75 TF	1	PTFE	1000 cm ² (155 in ²)	10-12 mm (3/8-1/2") SB	10-12 mm (3/8-1/2") SB	0.2	3	1
6710-7502	Polycap 75 TF	0.2	PTFE	1000 cm ² (155 in ²)	10-12 mm (3/8-1/2") SB	10-12 mm (3/8-1/2") SB	0.9	13	1
6711-7501	Polycap 75 TF	0.1	PTFE	1000 cm ² (155 in ²)	1/4" MNPT	6-10 mm (1/4-3/8") SB	1.7	25	1
6711-7502	Polycap 75 TF	0.2	PTFE	1000 cm ² (155 in ²)	1/4" MNPT	6-10 mm (1/4-3/8") SB	0.9	13	1

(1) SB: Stepped Hose Barb (dimensions refer to the I.D. of the hose)
MNPT: Male National Pipe Thread

(2) Typical Values
* Bubble Point and Flow Rate Values for 6711-7505 not available at time of printing

Legal.

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Whatman is a trademark of GE Healthcare companies.

© 2008–2013 General Electric Company – All rights reserved.
Previously published May 2008.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit

www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gelifesciences.com>

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

imagination at work