VWR[®] Microbiological Air Sampler

SAS Super DUO

SAS Super ISOLATOR

SAS Super DUO ISOLATOR

INSTRUCTION MANUAL

European catalogue numbers:

710-0867

710-0866

710-0943

710-0944

710-0945

710-0946

710-0977

710-0978

Version: 1.2

Issued: July 2016

Legal address of manufacturer

Europe

VWR International Srl

Via S. Giusto 85

20153 Milano

Italy

Tel.: + 39 02-3320311

<http://be.vwr.com>

Country of origin: Italy

Table of contents

Warning	5
Safety information.....	5
Unit descriptions.....	6
Technical features	7
Intended use	7
Principle	7
The basic idea.....	8
Brief instructions	8
Brief instructions for SAS Isolator	8
Brief operating instructions for DUO SAS Super 360 - DUO SAS Isolator	8
Practical use of contact plates	9
List of menus and utility sub menus	9
Preliminary inspection.....	10
Holder adjustment.....	10
Petri head adapter (optional)	10
90 mm Petri dish filling.....	10
Tripod installation (optional).....	10
Functions	10
Operations	10
Start with the same air volume as the previous sample	11
‘STANDARD MODE’	11
‘USER MODE’	12
‘PROGRAM MODE’	12
‘DELAY MODE’	13
‘MULTIMODE’	14
intervals).....	14
UTILITY MODE’	15
‘SET TIME’	16

'MODE'	16
'SET AUTOSWITCH'	16
'PRINT'	17
'DISPLAY RECORD'	17
'CLEAR RECORD'	18
'LANGUAGE'	18
'IDENTIFY'	18
'SAMPLING SITE'	18
'CALIBRATION DUE' (optional).....	19
List of system messages	19
Accessories	20
Troubleshooting	21
Technical service	22
Warranty	22
Compliance with local laws and regulations	22
Equipment disposal	23

Warning

This document is the property of VWR International S.r.l. - Milan - Italy

It may not be duplicated or distributed without the owner's authorisation.

Patent Pending

Manufactured by VWR International S.r.l. - Milan, Italy

Safety information

Use the device just for the indicated purposes.

The device must be correctly used according to this instruction manual, before starting any operation.

Immediately replace any electrical cable when damaged.

Always disconnect the charger before:

- Repairing or maintenance; these operations must be carried out by qualified staff
- Keep the unit clean

Use original spare parts and accessories for any replacement.

Do not use this device in the presence of explosive gas.

Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.



Be aware that the voltage and frequency of the electrical system are compatible with the electrical requirements of the battery charger.

Never use a non OEM charger to charge the air sampler.
Use of improper charger may damage the unit.

Unit descriptions

1. Command unit
2. Stainless steel head
3. Motor housing
4. "Clear" button
5. "Enter" button
6. Arrow "Up"
7. "Start" button

8. Arrow "Down"
9. Standing feet
10. Printer connection (not used)
11. Main switch
12. Serial number
13. CE mark
14. Charger connection

Technical features

Model	DUO SAS 360 Contact	DUO SAS 360 Petri	DUO SAS ISOLATOR Contact	DUO SAS ISOLATOR Petri
Cat. No.	710-0867	710-0866	710-0977	710-0978
Airflow	180 l/min	180 l/min	180 l/min	180 l/min
Used with	55 mm contact plates	90 mm Petri dishes	55 mm contact plates	90 mm Petri dishes
Portable	Yes	Yes	Yes	Yes
Battery	Rechargeable	Rechargeable	Rechargeable	Rechargeable
Battery life	40 000 l	40 000 l	40 000 l	40 000 l



Model	ISOLATOR 100 Contact	ISOLATOR 180 Contact	ISOLATOR 100 Petri	ISOLATOR 180 Petri
Cat. No.	710-0943	710-0944	710-0945	710-0946
Airflow	100 l/min	180 l/min	100 l/min	180 l/min
Used with	55 mm contact plates	55 mm contact plates	90 mm Petri dishes	90 mm Petri dishes
Portable	Yes	Yes	Yes	Yes
Battery	Rechargeable	Rechargeable	Rechargeable	Rechargeable
Battery life	70 000 l	40 000 l	70 000 l	40 000 l

Intended use

Principle

The Surface Air System (SAS) encompasses several models which use the same principle. Air is aspirated at a fixed speed for variable time through a head which has been machined with a series of small holes of a special design. The resulting laminar airflow is directed onto the agar surface of a “RODAC Plate” (or a Petri dish) containing medium for microbiological analysis.

When the selected sampling cycle is completed, the plate is removed and incubated. The organisms are then visible to the naked eye and can be counted for the assessment of contamination levels.

The basic idea

The key features of the Surface Air System (SAS) are:

1. To use a simple and inexpensive “contact plate” (“SURFAIR PLATE”, “RODAC”) for surface, hands or air control. These plates are very well known and easily available. They can be purchased ready poured with different media.
2. To sample a known volume of air for a variable time to provide a range of sampling volumes.
3. To aspirate air in a laminar flow pattern with sufficient velocity to impact organisms onto an agar surface.
4. To accumulate data on the level of hygiene in each environment so that fluctuations can be monitored.
5. To take advantage of advanced electronics for more reliable results in different operating conditions.
6. To have the flexibility to choose between 55 mm contact plates or 90 mm standard Petri dishes.
7. To organise sequential sampling to obtain a more representative sample under actual operating conditions.

Brief instructions

Brief instructions for SAS Isolator

The airflow is displayed after the automatic presentation every time the instrument is switched on.

Press ‘ON/OFF’ switch.

Press ‘START’ button to sample the same air volume as the last sampling cycle.

To change the air volume, use the arrow buttons and select “Standard Mode” from the menu.

Refer to the instruction manual to change settings.

Press ‘ENTER’ to confirm selection.

Press arrows for menu selection to modify other pre-set parameters:

STANDARD MODE
USER MODE
PROGRAM MODE
DELAY MODE
MULTI MODE
UTILITY MODE

Press ‘ENTER’ to reach the sub menu of the chosen parameter.

Refer to instruction manual to change settings.

Press ‘ON/OFF’ switch at the end of sampling operations.

Press ‘CLEAR’ each time you need to end an action. Then the unit will then return to the initial configuration.

Brief operating instructions for DUO SAS Super 360 - DUO SAS Isolator

Press the ‘ON/OFF’ switch (black button).

When the display will show the message ‘SELECT HEAD’, press ‘ENTER’.


Press the ‘UP’ or ‘Down’ arrows to select the ‘LEFT HEAD’ or the ‘RIGHT HEAD’ or ‘LEFT+RIGHT HEAD’.

Press ‘ENTER’ to confirm selection.

Press ‘START’ button to sample the same air volume as the last sampling time.

To change the volume of air or other parameters, follow the instructions for SAS ISOLATOR.

Practical use of contact plates

1. Remove the aspirating head.
2. Insert an identified, closed and prepared plate and remove the plate lid.
3. Replace the aspirating head.
4. Select required airflow and start the unit. The airflow is directed into the agar surface of the plate.
5. At the end of the cycle, remove the aspirating head.
6. Close and remove the plate.
7. Incubate.
8. Count the colonies, recording the results on the microbiological air sampling report and read the results.

List of menus and utility sub menus

START FOR XXX	This message indicates the air sampler is ready to sample the volume shown (XXX) after pressing the 'Start' button
USER MODE	With this function you can to modify the 8 programmable sampling volumes
UTILITY MODE	This leads to the following sub menu:
	SET TIME To adjust the date and time
	SAMPLING SITE To identify the sampling point
	IDENTIFY To identify the operator's name
	LANGUAGE To display text in different languages
	CLEAR RECORD To clear the recorded sampling data
MENU	DISPLAY RECORD To show the recorded sampling data

	AUTO SWITCH OFF To disconnect the automatic switch-off when the 'Infrared Remote' is used
	PRINT To export the recorded sampling data
	MODE To select the aspiration mode between single and multiple
	CALIBRATION DUE To check the calibration due date

Preliminary inspection

The apparatus is subject to specific quality tests before shipment and is carefully packed to avoid possible damage during transit. However, a visible check should be carried out as soon as possible to determine any damage. This must be reported immediately. The following procedures should be followed to check that the unit is working properly.

The battery pack of the instrument must be charged for at least 14 hours before running the test.

After the 'ON/OFF' button is pressed an automatic visual presentation will be displayed.

Holder adjustment

The contact plate and Petri dish holders can be adjusted (using a screw driver) in case the available plates is slightly different in diameter from the standard 55 mm contact plate or 90 mm Petri dish.

Petri head adapter (optional)

An adapter is available in aluminium or stainless steel which allows the use of standard 90 mm Petri dishes with SAS for contact plates.

90 mm Petri dish filling

We recommend you fill the standard 90 mm disposable Petri dish with no more than 18 - 20 ml of agar to avoid the medium touching the inside surface of the aspirating head.

Tripod installation (optional)

The SAS can be fixed to a table, floor tripod or wall support. The screw thread connection is located under the unit, between the two front feet.

When not in use, the connection cavity is protected by a metal screw insert.

Functions

Operations

Press the 'ON/OFF' switch (black button) to switch on the unit.

The 'ON/OFF' switch (black button) must be pushed twice to switch the unit on if the automatic switch-off timer was used previously.

Press 'CLEAR' at any time to return to the initial 'START FOR' main screen.

The air sampler is provided with eight fixed air volumes ('Standard Mode') stored in the memory and eight programmable air volumes ('USER MODE') selectable by the operator.

Volume selection is made by pressing the "▲" or "▼" arrows when the programme is in the relevant sub menu.

Suggested air volumes:

- Contaminated areas (communities, processing rooms, etc) 10 - 200 litres of air
- Normal areas (laboratory benches, houses, etc) 200 - 500 litres of air
- Sterile or high risk areas (cleanrooms, operating theatres, etc) 500 - 1000 litres of air.

Start with the same air volume as the previous sample

SAS SUPER ISOLATOR

Switch on the SAS using the 'ON/OFF' switch. You will then see the following display with the last aspirated volume of air:

```
START FOR   XXX
MM/DD/YY - HH:MM
```

Push the 'START' button to start a sampling cycle with the displayed volume of air.

```
[ XX]      →   XXX
MM/DD/YY - HH:MM
```

DUO SAS SUPER 360 - DUO SAS Isolator

Switch on the SAS using the 'ON/OFF' switch. You will then see the following display with the last aspirated volume of air:

```
SELECT HEAD
```

Select the desired head by using the "▲" or "▼" arrows and press 'ENTER', the main screen will be displayed

```
START FOR   XXX
MM/DD/YY - HH:MM
```

Push the 'START' button to start a sampling cycle with the displayed volume of air.

```
[ XX]      →   XXX
MM/DD/YY - HH:MM
```

'STANDARD MODE'

The standard mode function allows the operator to select a sampling volume from a list of eight standard volumes.

Switch on the instrument and wait until the main screen is displayed:

```
START FOR   XXX
MM/DD/YY - HH:MM
```

Press the "▼" or "▲" arrows to reach the

```
STANDARD MODE
MM/DD/YY - HH:MM
```

Press the 'ENTER' button, now you can scroll the eight sampling volume. The selectable volumes are 10, 20, 30, 50, 100, 200, 500, 1000 litres).

S.PROG XXXX
MM/DD/YY - HH:MM

Once the desired sampling volume has been reached press the 'ENTER' button to confirm. The instrument is now ready to sample the selected volume of air. Press the 'START' button to begin the sampling cycle

START FOR XXX
MM/DD/YY - HH:MM

'USER MODE'

If the operator wish to use volumes other than the standards volumes, up to 8 additional volumes can be programmed into the unit.

The user-selectable volumes can be selected from the USER MODE, while the selectable volumes can be modified using the PROGRAM MODE function describe in the next paragraph. Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the "▼" or "▲" arrows to reach the

USER MODE
MM/DD/YY - HH:MM

Press the 'ENTER' button, now you can scroll the eight user-defined sampling volumes.

S.PROG XXXX
MM/DD/YY - HH:MM

Once the desired sampling volume has been reached press the 'ENTER' button to confirm. The instrument is now ready to sample the selected volume of air. Press the 'START' button to begin the sampling cycle

START FOR XXX
MM/DD/YY - HH:MM

The total number of available 'USER MODE' programmes is eight. The maximum volume of air for each sampling cycle is 1999 litres.

'PROGRAM MODE'

Using this procedure, it is possible to memorise up to eight different volumes (from 1 to 1999 litres of air). (See chapter 4.4.).

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the “▼” or “▲” arrows to reach the

PROGRAM MODE
MM/DD/YY - HH:MM

Press the ‘ENTER’ button, now you can scroll the eight user-selectable sampling volumes.

S.PROG XXXX
MM/DD/YY - HH:MM

Once the sampling volume to be modified has been reached, press the ‘ENTER’ button to confirm.

The display will show the old and the new volume.

XXX → 0000
MM/DD/YY - HH:MM

Press the “▼” or “▲” arrows to select the first digit and press ‘ENTER’ to confirm

XXX → X000
MM/DD/YY - HH:MM

Press the “▼” or “▲” arrows to select the second digit and press ‘ENTER’ to confirm. Repeat the previous step until the last digit is reached. When the last digit is confirmed, the new volume will be memorised and added to the list of eight user-selectable volumes. If you want to start a sampling cycle with the new value you have to enter the USER MODE and select it by using the procedure described in the USER MODE paragraph.

‘DELAY MODE’

With this function it is possible to delay the beginning of the sampling cycle. If a delay has been selected, the sampling cycle will begin after the selected delay time.

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the “▼” or “▲” arrows to reach the

DELAY MODE
MM/DD/YY - HH:MM

Press the ‘ENTER’ button, the delay time is displayed

DELAY 01 MIN
MM/DD/YY - HH:MM

Select the desired delay time with the “▼” or “▲” arrows. The selectable delay values are 1, 2, 3, 5, 10 or

20 minutes. Press 'ENTER' button to confirm.

When a delay has been selected, after the 'START' button has been pressed, a delay warning message will be displayed

** D E L A Y **

'MULTIMODE'

This programme is very useful for extending the time of sampling with the purpose of obtaining a more representative environmental sample "in actual operating conditions". The total air volume to be sampled is aspirated with two or more sub volume aspirations (e.g.:1000 litres in ten runs of 100 litres at five minute intervals).

Before entering 'MULTI MODE' you should decide:

- (a) Total volume of air to be sampled onto the contact plate.
- (b) Number of runs.
- (c) Interval time between runs.

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the "▼" or "▲" arrows to reach the

MULTI MODE
MM/DD/YY - HH:MM

Press the 'ENTER' button, the interval time is displayed

INTERVAL TIME
 05 MIN

Select the desired interval time between runs with the "▼" or "▲" arrows. The selectable interval values are 5, 10, 15, 20, 25, 30, or 60 minutes.

Press the 'ENTER' button to confirm the interval time, the numbers of runs is then displayed

NUMBER OF RUNS
 20 NC

Select the desired number of cycles with the "▼" or "▲" arrows. The selectable values are 2, 3, 4, 5, 6, 7, 8, 9, 10, 15 or 20 runs.

Press the 'ENTER' button to confirm the number of runs, the single run volume is then displayed.

SINGLE RUN VOL.
XXXX

Select the desired digit using the "▼" or "▲" arrows then press the 'ENTER' button to skip to next digit.

When all the digits have been selected the total volume is displayed.

TOT. L. XXXX

Press the 'ENTER' button to confirm.

NOTE: Before running multimode sampling, set the aspirating mode to MULTIMODE (see UTILITY MODE paragraph)

Starting a multimode sampling cycle

Before running multimode sampling, set the aspirating mode to MULTIMODE (see chapter 4.8.2).

When the sampler aspiration mode is set to MULTIMODE, the main display will alternately show the number of cycles and the volume of each cycle:

START FOR
MM/DD/YY - HH:MM

NC 02x1000 L.
MM/DD/YY - HH:MM

Press the 'START' button to run multiple sampling, the display will show the aspirating volume in the first row and the residual number of the cycle in the second row:

[XX] → XXX
- 02 CYCLES

At the end of the run, during the interval, the number of "residual cycles" (runs) that are remaining are shown on the display:

- 02 CYCLES
** D E L A Y **

UTILITY MODE'

With this function it is possible to modify all the settings of the air sampler. The UTILITY MODE is divided into ten sub menus:

Set time; Mode; Set Autoswitch; Print; Display Record; Clear Record; Language; Identify; Sampling site; Calibration due.

Always start from UTILITY MODE to reach one of these sub menus. Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the "▼" or "▲" arrows to reach the

UTILITY MODE
MM/DD/YY - HH:MM

Press the 'ENTER' button, the first sub menu is listed.

‘SET TIME’

From the UTILITY MODE select the ‘SET TIME’ function

This option is used to programme day, month, year and time of the day.

SET TIME

MM/DD/YY - HH:MM

Press ‘ENTER’ and “▼” or “▲” arrows to change month

Press ‘ENTER’ and “▼” or “▲” arrows to change day

Press ‘ENTER’ and “▼” or “▲” arrows to change year

Press ‘ENTER’ and “▼” or “▲” arrows to change hours

Press ‘ENTER’ and “▼” or “▲” arrows to change minutes

Press ‘ENTER’ to confirm and exit.

‘MODE’

From the UTILITY MODE select the ‘MODE’.

This option is used to select the aspirating mode between SINGLE MODE and MULTIMODE.

MODE

MM/DD/YY - HH:MM

Press ‘ENTER’ and “▼” or “▲” arrows to change between SINGLE MODE and MULTIMODE

MULTI MODE

MM/DD/YY - HH:MM

SINGLE MODE

MM/DD/YY - HH:MM

Select the desired aspirating mode and press ‘ENTER’ to confirm and exit.

This setting will be stored in memory even after the instrument is switched off.

‘SET AUTOSWITCH’

To save battery consumption the SAS automatically switches off after 4 minutes. The SET AUTOSWITCH option is used to enable or disable the automatic switch-off.

From the UTILITY MODE select the ‘SET AUTOSWITCH’

SET AUTOSWITCH

MM/DD/YY - HH:MM

Press ‘ENTER’ and “▼” or “▲” arrows to change between the two options.

AUTOSWITCH ON

MM/DD/YY - HH:MM

AUTOSWITCH OFF
MM/DD/YY - HH:MM

Select the desired options and press 'ENTER' to confirm and exit.

This setting will be lost after the instrument is switched off.

'PRINT'

The SAS can be connected with a RS232 cable to the SAS printer (Cat. No. ZZ30PCDPT100-I-12) or to a Personal Computer with the SAS software (Cat. No. 710-0975). This function is used to start the data transfer from the SAS to an external device.

From UTILITY MODE select the 'PRINT'

PRINT
MM/DD/YY - HH:MM

Press 'ENTER' to start sending the data, the display will show:

SENDING DATA

For more information see the printer or software user manual.

'DISPLAY RECORD'

The last 99 samples are memorised in the file 'DISPLAY RECORD'. Each sample is identified in chronological order and shows the date, time, operator, site and volume of air sampled.

From the UTILITY MODE select the 'DISPLAY RECORD'.

DISPLAY RECORD
MM/DD/YY - HH:MM

Press 'ENTER' to see the memorised values

001 MM/DD HH:MM
ID. /SITE VOL. OP

The following parameters are recorded:

Progressive number; Month and Day; Time; Identification / Site; Sampling Volume; Option

If there is no data stored in memory the display will show

* * * * *
* * * * *

‘CLEAR RECORD’

This option is used to delete all the data memorised in the DISPLAY RECORD. Before starting this procedure, please be certain that existing data is not required or that it has been downloaded.

From the UTILITY MODE select the ‘CLEAR RECORD’.

CLEAR RECORD

MM/DD/YY - HH:MM

Press ‘ENTER’ to delete all the stored data, the display will show the resetting status

RESETTING...

▶▶▶▶

‘LANGUAGE’

Menu text can be selected choosing from different languages. From the UTILITY MODE select the ‘LANGUAGE’.

LANGUAGE

MM/DD/YY - HH:MM

Press ‘ENTER’ to list all the available languages,

ENGLISH

Use the “▼” or “▲” arrows to select the desired language and press ‘ENTER’ to confirm.

There are six available languages:

English, French, Spanish, German, Portuguese and Italian.

‘IDENTIFY’

This option is used to identify the operator. This should be changed if different operators use the sampler and especially if the data is to be printed.

From the UTILITY MODE select the ‘IDENTIFY’.

IDENTIFY

MM/DD/YY - HH:MM

Press ‘ENTER’ to modify the identification code.

IDENTIFY

XXXX

Select the desired character using the “▲” or “▼” arrows and confirm with ‘ENTER’.

Select the next desired characters and confirm them by pressing ‘ENTER’.

The selected values are stored in the memory and it will be kept also after the instrument has been switched off.

‘SAMPLING SITE’

The site identification should be changed for samples taken at different sites especially if the results are to be printed.

From the UTILITY MODE select 'SAMPLING SITE' function.

SAMPLING SITE
MM/DD/YY - HH:MM

Press 'ENTER' to modify the sampling site

SAMPLING SITE
XXXX

Select the desired character using the "▲" or "▼" arrows and confirm with 'ENTER'.

Select the next desired characters and confirm them by pressing 'ENTER'.

The selected values are stored in the memory and it will be kept also after the instrument has been switched off.

'CALIBRATION DUE' (optional)

The SAS air sampler can be equipped with a calibration reminder feature. With this function the user can see when the next calibration should be performed.

From the UTILITY MODE select the 'CALIBRATION DUE'.

CALIBRATION DUE
MM/DD/YY - HH:MM

Press 'ENTER' to see the next calibration date:

NEXT CAL
MM/YYYY

The display shows the month in which the calibration should be performed.

One month before the calibration is due; the SAS will display the following warning message:

CALIBRATION DUE
MM/YYYY

When this message is displayed, press 'ENTER' to acknowledge.

When the calibration has expired, the SAS will display the following warning message:

CAL EXPIRED

When this message is displayed, press 'ENTER' to acknowledge.

The calibration timer can be reset only by a qualified technician.

List of system messages

Besides the messages described in the previous paragraph, the following messages can be displayed:

LOW BATTERY	Battery must be recharged
MOTOR ERROR	Motor not working or not connected

- OFF	The instrument is switching off
NEXT CAL	The calibration due date is approaching
CAL EXPIRED	The calibration due date has expired

Accessories

Description	Cat. No.
Aspirating head for contact plates, Ø 55 mm	
Stainless steel aspirating head for contact plates, Ø 55 mm	710-0880
Sterile daily head for contact plates, Ø 55 mm	710-0890
Aluminium aspirating head for contact plates, Ø 55 mm	710-0892
Aspirating head for Petri dishes, Ø 90 mm	
Stainless steel aspirating head for Petri dishes, Ø 90 mm	710-0878
Sterile daily head for Petri dishes, Ø 90 mm	710-0891
Aluminium aspirating head for Petri dishes, Ø 90 mm	710-0886
Accessories	
SAS Super universal battery charger	710-0993
Remote control	710-0962
Soft carrying case for DUO SAS	710-0888
Aluminium carrying case for DUO SAS	710-0876
Aluminium carrying case for SAS ISOLATOR	710-1022
Floor tripod	710-0889
SAS-Holder table and wall stainless steel	710-0963
Adaptor* to convert contact plate model to accept 90 mm Petri dishes	710-0882
SAS stainless steel Petri head + adaptor	710-0877
SAS aluminium Petri head + adaptor	710-0879

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
Unit does not start	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and replace if required (Cat. No. 710-0993)
	Battery is too old	Check battery and replace if required
Battery discharges after few minutes working	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and replace if required (Cat. No. 710-0993)
	Battery is too old	Check battery and replace if required
“LOW BATTERY” message	Battery is low	Recharge battery
“CAL EXPIRED” message	Instrument needs to be calibrated	Send the instrument to VWR International VWR or an authorised dealer
Infrared remote switch doesn't switch on the unit	Remote control battery is low	Replace infrared remote battery
Microbiological media is dehydrated after sampling	Media is corrupted	Check expiration date of media and if the agar was not dehydrated before sampling
	Sampling time is too long	Shorten the sampling time

Technical service

Web resources

Visit the VWR website at vwr.com for:

- Complete technical service contact information
- Access to the VWR online catalogue and information about accessories and related products
- Additional product information and special offers

Contact us for information or technical assistance contact your local VWR representative or visit vwr.com

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handing it over to an authorised facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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VWR International AB
VWR International AG
Teknolojileri Ltd.Şti.
VWR International Ltd