

VWR® GEL DOCUMENTATION SYSTEMS

VWR® Basic

VWR® Smart3/Smart3 EZ

VWR® Imager2

VWR® Imager CHEMI
Premium

VWR® CHEMI *only*



Choose from the cost effective VWR® Basic system to the most advanced VWR® Imager CHEMI Premium for fluorescence and chemiluminescent imaging. The range includes models adapted to many different applications so you can choose the best solution for your work.

View the potential configurations below.

	Basic	Smart3/Smart3 EZ	CHEMI only	Imager2	Imager CHEMI Premium
UV light transillumination	✓	✓		✓	✓
Blue light transillumination	✓	✓		✓	✓
R/G/B and IR epi illumination					✓
UV light epi illumination				✓	✓
White light epi illumination		✓	✓	✓	✓
Converters (UV to white or blue light)	✓	✓		✓	✓
Colorimetric gels (e.g. Coomassie, silver stains)	✓	✓		✓	✓
Colony plates	✓	✓		✓	✓
Slot blots			✓		✓
Chemiluminescence			✓		✓
Bioluminescence			✓		✓
Stain-free					✓



VWR® Basic

Cost effective, quality imaging

This entry level system can be used with UV, blue or visible light transilluminators. These flexible lighting options make the VWR® Basic system suitable for imaging a wide range of samples.

- High resolution digital camera for high quality images with precise band separation
- Full PC interface enables control of camera with real time imaging
- Programmable capture settings

Using the VWR Basic gel documentation system you can produce images of gels stained with many fluorescent and colorimetric dyes, including:

- Coomassie Blue
- Silver stain
- Ethidium bromide
- GelRed™
- GelGreen™
- EZ-VISION
- SYBR® Gold/Green/Safe
- SYPRO® Red/Ruby
- Fluorescein

As new dyes are released, we work to optimise their use with the Basic system. Contact our technical support team if you need information on imaging dyes.

VWR® Basic	Technical specifications
Camera	
Resolution	18 MP
Sensor/image depth	24 bit (colour)
Lens	Manual zoom, 3x optical
Power	Rechargeable battery (supplied)
Data type	JPG
Interface	USB 2.0
Image storage	Mass storage
Viewing area	20x20 cm
Illumination	
Compatible with VWR stand-alone UV, blue or white light transilluminators with max. 20x20 cm filter size	Sold separately
Darkroom	
Hood	Sits directly on transilluminator
Stand-alone and PC-operated (optional)	PC sold separately
Image viewing area	20x20 cm
Filter(s)	
UV/blue filter	Included
Software	
VWR® Image Capture Software	Included
VWR® Gel Documentation Software (for analysis)	Included (1 license)
Dimensions	
HxDxW	37x21x28 cm

Description	Includes	Cat. No.
VWR® Basic	18 MP camera, hood, UV filter, 1 copy of the VWR Gel Documentation Software. Transilluminator and converter not included	730-1470

VWR® Smart3/Smart3 EZ

Quick and easy image capture

The VWR® Smart3/Smart3 EZ gel documentation systems come with an integrated PC and large touch display. Icon-driven menus guide you through the various functions to capture, print and save images with one simple touch on the large screen.

The sliding transilluminator offers easy access to your sample. The filter drawer, with interchangeable filters (optional), enables detection of a wide range of different fluorescent stains. The UV to white light converter screen enables imaging of conventional and colorimetric samples, including Coomassie gels, silver stains, autoradiographs and colony plates.

The camera has an excellent resolution of 3 MP and images can be rapidly saved as TIFF, BMP or JPEG files onto a USB drive or directly to a PC via a network. The VWR® Smart3 system has a manual zoom lens and the Smart3 EZ is equipped with a motor-driven zoom lens.

Smart3 systems include the innovative VWR Gel Documentation Software that makes it easy to analyse and quantify your samples. Intuitive handling and option to automate analysis, save time and reduce workload.



VWR® Smart3/Smart3 EZ gel documentation system can produce images of gels stained with a wide range of fluorescent and colorimetric dyes, including:

- Coomassie Blue
- Silver stain
- Ethidium bromide
- GelRed™
- GelGreen™
- EZ-Vision®
- peqGREEN
- SYBR® Gold/Green/Safe
- SYPRO® Red/Ruby
- Fluorescein
- Pro-Q® Diamond
- Deep Purple™

As new dyes are released, we work to optimise their use with the VWR Smart systems. Contact our technical support team if you need information on imaging dyes.

VWR® Smart3/Smart3 EZ	Technical specifications
Camera	
Resolution	3 MP
Sensor/image depth	12/16 bit (4,096/65,536 greyscale)
Lens	Manual zoom, f/1,4 (Smart3) Motor-driven zoom, f/1,4 (Smart3 EZ)
Illumination	
VWR UV transilluminators, 20x20cm	Sold separately
UltraSlim LED blue light transilluminator, 10x12cm	Sold separately
Epi white light	Yes
Visible light converter (UV to white light)	Sold separately
Blue light converter screen (UV to blue light)	Sold separately
Darkroom	
Cabinet	Stand-alone with integrated PC and display (no external PC needed)
Image viewing area	20x20 cm
Filter(s)	
UV/blue filter	Included
Software	
VWR® Image Capture Software	Included
VWR® Gel Documentation Software (for analysis)	Included (2 licenses)
Dimensions	
HxDxW	51x39x46,5 cm

Description	Cat. No.
VWR® Smart3 Complete system, 3 MP 12/16 bit camera, manual zoom lens, 2 copies VWR® Gel Documentation Software, darkroom, network connection, built-in processor and touch screen (no external PC required). Transilluminator and converter not supplied.	730-1459
VWR® Smart3 EZ Complete system, 3 MP 12/16 bit camera, motor-driven zoom lens, 2 copies VWR® Gel Documentation Software, darkroom, network connection, built-in PC and touch screen (no external PC required). Transilluminator and converter not supplied.	730-1460

VWR® CHEMI *only* can produce images of Western blots stained with chemiluminescence substrates, including:

VisiGlo™ HRP chemiluminescent substrates
ECL Western blotting substrates

VWR CHEMI <i>only</i>	Technical specifications
Camera	
Image resolution	4 MP
Sensor/image depth	16 bit (65,536 greyscale)
Lens	Fixed, f/0,95
QE @ 425 nm	0,73
Cooling	DT -57 °C
Illumination	
Epi white light	Yes
Use with external PC	Optional
Darkroom	
Cabinet	PC-operated (PC sold separately)
Image viewing area	11x8 cm
Software	
VWR Image Capture Software	Yes
VWR Gel Documentation Software	Yes (2 licenses)
Dimensions	
HxDxW	37,4x47x46 cm

VWR® CHEMI *only*

Chemiluminescence imaging

The new VWR® CHEMI *only* is dedicated to chemiluminescence imaging. This system is built for high performance and automation featuring a next generation high quantum efficiency CCD camera for even greater sensitivity. VWR® CHEMI *only* allows automatic capture of a quality image of any Western blot with a single 'click'.

Western blot imaging has never been as easy as with the VWR® CHEMI *only*. The system automatically selects the right imaging conditions for any blot regardless of the chemiluminescent reagents being used. All chemiluminescence applications can easily be handled by the VWR® CHEMI *only*. Produce superb images even from the faintest of signals.

The compact instrument sits on any laboratory bench taking up very little space.

The VWR® CHEMI *only* is ideal for small blots (up to 11x8 cm). The cooled 16 bit camera (ΔT -57 °C) with fixed lens (f/0,95) and data feedback amazes with a great quantum efficiency (73% @ 425 nm) and high resolution (4 MP).



Description	Includes	Cat. No.
VWR® CHEMI <i>only</i>	4 MP image resolution, 12/16 bit, 2 copies of VWR® Gel Documentation Software, darkroom, 730-1471 Image Capture Software.	730-1471

VWR® Imager systems

Fluorescence and chemiluminescence imaging

VWR Imager systems offer a wide range of different options depending on the configuration. Suitable for colorimetric, fluorescence, chemiluminescence and bioluminescence applications.

Choose from the following models:

- VWR® Imager2 for colorimetric and fluorescence applications only
- VWR® Imager CHEMI Premium for standard applications plus chemiluminescence and bioluminescence

Each system features a CCD camera, a motorised 7-position filter wheel and a motor-driven zoom lens to allow for a wide range of imaging applications.

The Imager CHEMI Premium comes with a cooled camera (ΔT -57 °C) that also enables use in chemiluminescence and bioluminescence applications.

Choose from a wide range of lighting options for both transillumination and epi illumination including UV, blue and white light.

The VWR® Imager's simple design allows for a wide range of configurations. Select from two cameras, many different filters and lighting options. The darkroom is fully microprocessor controlled, so that functions such as camera settings, lens control (with data feedback on VWR® Imager CHEMI Premium), filter selection and lighting can all be set up via PC (not included). In addition, optional settings such as manual or auto exposure time, neutral fielding and extended dynamic range are easily accessible for the user. All settings can be saved in user-defined configurations.

Technical specifications	VWR® Imager2	VWR® Imager CHEMI Premium
Camera		
Image resolution	3,8 MP	4 MP
Sensor/image depth	12/16 bit (4,096/65,536 greyscales)	16 bit (65,536 greyscales)
Lens	Motor-driven zoom, f/1,2	Motor-driven zoom, f/1,2 with data feedback
QE @ 425 nm	N/A	0,73
Cooling	No	DT -57 °C
Illumination		
VWR UV and blue light transilluminators	Sold separately	Sold separately
Epi white light	Yes	Yes



VWR® Imager2 gel documentation systems can produce images of gels stained with many fluorescent and colorimetric dyes, including:

- Coomassie Blue
- Silver stain
- Ethidium bromide
- GelRed™
- GelGreen™
- EZ-VISION
- peqGREEN
- SYBR® Gold/Green/Safe
- GelStar®
- SYPRO® Red/Ruby/Orange
- Fluorescein
- Rhodamine Red™
- Texas Red®
- Pro-Q® Diamond
- Deep Purple™
- GFP

VWR® Imager CHEMI Premium also offers imaging of:

- Western blots
- Chemiluminescence
- Bioluminescence
- Stain-free gels
- Alexa Fluor® dyes
- CF™ dyes
- Cy® dyes

We work with many more dyes than listed above, and as new dyes are released, we work to optimise their use with the VWR Imager range so please ask your VWR contact for updates.

Computer (not supplied)

The darkroom is fully microprocessor controlled. Camera settings, lens control, filter selection and lighting are set-up via a PC. User-configurations can then be saved for future use. Please ask us for minimum computer requirements.

Technical specifications	VWR® Imager2	VWR® Imager CHEMI Premium
Epi UV light modules	Sold separately	Sold separately
Epi R/G/B and IR LED modules	No	Sold separately
Visible light converter (UV to white light)	Sold separately	Sold separately
Blue light converter screen (UV to blue light)	Sold separately	Sold separately
Darkroom		
Cabinet	PC-operated (PC sold separately)	PC-operated (PC sold separately)
Filter wheel	Yes (7-position, motor-driven)	Yes (7-position, motor-driven)
Image viewing area	25,5x21 cm	30,5x22,7 cm
Software		
VWR® Image Capture Software	Yes	Yes
VWR® Gel Documentation Software (for analysis)	Yes (2 licenses)	Yes (2 licenses)
Dimensions		
HxDxW	84x45x57 cm	84x45x57 cm

Description	Cat. No.
VWR® Imager2 3,8 MP 12/16 bit camera, 7-position filter wheel, 2 copies of VWR® Gel Documentation Software, darkroom, VWR® Image Capture Software. Transilluminator and converter not supplied.	730-1458
VWR® Imager CHEMI Premium 4 MP image resolution, 12/16 bit, 7-position filter wheel, 2 copies of VWR® Gel Documentation Software, darkroom, LED gantry, VWR® Image Capture Software. Transilluminator and converter not supplied.	730-1469

Cameras

Both VWR® Imager systems are fitted with CCD cameras. The VWR® Imager2 is equipped with a 3,8 MP camera, enabling analysis of e.g. colorimetric and fluorescent gels. The VWR® CHEMI Premium offers a 16-bit cooled (ΔT -57 °C) CCD camera with a resolution of 4 MP, allowing imaging of colorimetric, fluorescence, chemiluminescence and bioluminescence applications.

Lenses

All VWR® Imager systems record lens settings in order to meet GLP requirements. They are equipped with fast f/1,2 zoom lenses and directly controlled by the VWR® Image Capture Software. Both VWR Imagers deliver a superb optical performance.

Darkroom

The units incorporate a fully light-tight darkroom best suited for chemiluminescent and fluorescent imaging (e.g. UV/blue transillumination and UV, R/G/B and IR epi illumination). With a wide opening door and a sliding transilluminator (optionally UV, blue or white light), VWR Imagers offer maximum convenience. The electronic door locks automatically if exposures over 30 seconds are selected to prevent accidental interruption of imaging proceedings. VWR Imager units are built in a modular form to allow a wide choice of accessories to extend applications. The VWR® Imager2 is also fully upgradable for chemiluminescence applications.

Lighting

Range of lighting options

- UV transilluminators (254, 302 or 365 nm or a combination of any 2 wavelengths, either with 20x20 cm or 25x30 cm illumination area)
- Blue light transilluminator for safe dyes such as GelRed™, GelGreen™, SYBR® dyes and EZ-Vision (20x16 cm illumination area)
- Visible light converter screen enabling white light transillumination for colorimetric applications such as Coomassie Blue and silver stains
- Blue light converter screen enabling blue light transillumination for 'Safe' dyes such as GelRed™, GelGreen™, SYBR® dyes and EZ-Vision. UV light epi illumination (254, 302 and 365 nm)
- White light epi illumination (standard for all VWR Imagers)
- Epi illumination for R/G/B and IR fluorescence multiplex-imaging (only for VWR® Imager CHEMI Premium. Fitting LED modules can be added at any time in the future)
- Neutral fielding screens

Filters

VWR Imager systems are fitted with an orange UV filter (572 – 625 nm) as standard. There are a number of additional filters that cover an extensive list of applications.

- Short pass filter - 516 – 600 nm
- Long pass filter - 611 – 641 nm
- Far right red filter - 670 – 780 nm
- Filters for R/G/B multiplex imaging
- Filters for Q-Dots

For more information please contact your local VWR sales office.




VWR® Image Capture Software

An advanced image capture software package specifically developed to simplify the process of capturing images. Every VWR Imager system includes VWR® Image Capture Software.

This fully automatic package controls camera integration, exposure time, lens and capture options, as well as auto-focus configurations. A complete database of dyes and reagents at your fingertips. A comprehensive post-processing tool box allows the user to add text, lines, circles, shapes and arrows to the captured images. One button technology – a single button is all that is required to capture the perfect image every time. Lens control for operating the aperture, focus and zoom position for the motor-driven lens (with displayed feedback data).

- Full on-screen control of darkroom, camera and lens functions
- Auto switch-off function for UV lighting
- 'Live image' feature in real time for easy focusing and positioning
- Automatic exposure control even for chemiluminescent imaging
- Series capture' for a series of images over different time periods (e.g. for determination of ideal exposure times for chemiluminescent samples)
- Saturation detection for quantitative applications
- Image annotations, enhancement and manipulation
- Image sharpening and smoothing
- Style control (lines, size, colour, fonts)
- Image cropping
- Save images in SGD private file format for true GLP
- Export images in various formats (TIF, BMP, JPEG)
- Direct link to the VWR® Gel Documentation Software for analysis



VWR® Gel Documentation Software

VWR® Gel Documentation Software is an advanced analysis software for use with any VWR Imager system. This highly automated software can rapidly analyse an image of the gel in a matter of seconds. Simplicity is the key feature of the VWR® Gel Documentation Software, results can be obtained in a click of a button and with minimal training.

Features of the VWR® Gel Documentation Software

ID analysis

- MW/BP analysis
- Automatic lane and peak detection
- Single or multiple gel analysis
- Automatic compensation for smiling or distorted bands – no skewing required
- MW/BP/Quantity calibrations
- Extensive reporting
- Spot analysis
- Spot grid – automatic
- Spot thresholding
- Quantitative and incidence analysis
- Multiple background subtraction methods

Colony counting

- Advanced colony counting algorithm
- Automatic counting/manual detection
- Two colour counting
- Automatic separator
- Light/dark colony selection
- Exclude regions function
- Sensitivity selector
- RGB multiplex analysis




Extensive applications and dye imaging

- *RGB multiplex analysis
- Quantitative dots (Q-dots)
- SYBR® dyes
- Cy® dyes
- Alexa Fluor® dyes
- SYPRO® dyes
- Pro-Q® dyes
- Plus many more

Visible light (for colorimetric applications such as Coomassie Blue, silver staining and various plate-based assay)

For applications such as Coomassie and plates.

- Optional conversion screen





730-1461
Transilluminator 20x20 cm, 302 nm



730-1389
Transilluminator Blue, 20x16 cm, 470 nm



730-1390
Transilluminator Blue Slim 10x12 cm, 470 nm



730-1260
Thermal printer, digital

VWR® transilluminators

All UV transilluminators slide out for easy access and offer homogeneous illumination for the best performance.

Main features and choices

- Safety cut-off (when in darkroom)
- Slide out of darkroom for easy access and sample positioning
- Different single and dual wavelengths available
- Different filter sizes available

Ordering information

Transilluminators for VWR® Basic

Transilluminator Blue slim is small and less expensive. 10x12cm sample size. Transilluminator Blue is larger and works best with the VWR® Imager2 and VWR® Chemi Premium darkrooms. Plus can be used with the VWR® Basic as this is bench top unit

Description	Cat. No.
Transilluminator 20x20 cm, 302 nm	730-1464
Transilluminator 20x20 cm, 365 nm	730-1465
Transilluminator 20x20 cm, 254 nm	730-1466
Transilluminator Blue, 20x16 cm, 488 nm	730-1467
Visible light converter (UV to white light)	730-1395

VWR transilluminators for VWR Smart3/Smart3 EZ

Description	Cat. No.
Transilluminator 20x20 cm, 302 nm	730-1461
Transilluminator 20x20 cm, 365 nm	730-1462
Transilluminator 20x20 cm, 254 nm	730-1463
Transilluminator Blue slim, 10x12 cm, 470 nm	730-1468
Visible light converter (UV to white light)	730-1395
Blue light converter (UV to blue light)	730-1493

Transilluminators for Imager2 and CHEMI Premium

Description	Cat. No.
Transilluminator 20x20 cm, 302 nm	730-1482
Transilluminator 25x30 cm, 302 nm	730-1483
Transilluminator 20x20 cm, 365 nm	730-1472
Transilluminator 25x30 cm, 365 nm	730-1473
Transilluminator 20x20 cm, 302/365 nm	730-1474
Transilluminator 25x30 cm, 302/365 nm	730-1475
Transilluminator 20x20 cm, 254 nm	730-1476
Transilluminator 25x30 cm, 254 nm	730-1477
Transilluminator 20x20 cm, 254/365 nm	730-1478
Transilluminator 25x30 cm, 254/365 nm	730-1479
Transilluminator 20x20 cm, 254/302 nm	730-1480
Transilluminator 25x30 cm, 254/302 nm	730-1481
Transilluminator Blue, 20x16 cm, 470 nm	730-1389
Transilluminator Blue Slim, 10x12 cm, 470 nm	730-1390
GX-Convert5 visible light converter (UV to white light)	730-1396
Blue light converter 25x30cm	730-1494
Blue light converter 21x26cm	730-1493
White light pad with brackets, 20x14 cm	730-1391

Description	Cat. No.
Epi UV module 254 nm	730-1394
Epi UV module 302 nm	730-1393
Epi UV module 365 nm	730-1392
Only for VWR® Imager CHEMI Premium	
Epi-LED blue module	733-2369
Epi-LED blue multiplexing module	733-2314
Epi-LED red module	733-2371
Epi-LED red multiplexing module	733-2372
Epi-LED green module	733-2370
Epi-LED green multiplexing module	733-2315
Epi-LED IR 740 multiplexing module	733-2316
Filter for epi-LED module, 800 nm (range 809 - 876; LiCor IRDye800), Chemi systems only	730-1524
525 nm filter for use with blue multiplexing LEDs	733-2305
Filter for Licor multiplexing 800 nm (809 - 876 nm)	730-1524
605 nm filter for use with green multiplexing LEDs	733-2480
705 nm filter for use with red multiplexing LEDs	733-2481

Printer options

VWR® Smart3, VWR® Smart3 EZ, VWR® Imager2, VWR® Imager CHEMI Premium, VWR® CHEMI *only* can all be used with the following digital thermal printer and paper.

Description	Cat. No.
Thermal printer, digital	730-1260
Thermal paper, matt	730-2892
Thermal paper, glossy	733-2000

VWRbioMarke™

For Life Scientists

OUR EXCLUSIVE PROGRAMME FOR YOUR SUCCESS IN PROTEOMICS - GENOMICS - CELL BIOLOGY

Instrumentation for many techniques
Consumables for cell culture, molecular biology and immunology
Research reagents, biological buffers, standards, media and stains
For more details contact your local VWR sales office.

Austria

VWR International GmbH
Graumannsgasse 7
1150 Vienna
Tel.: +43 1 97 002 0
Fax: +43 1 97 002 600
Email: info.at@vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: 016 385 011
Fax: 016 385 385
Email: vwr.be@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
Fax: +420 321 570 320
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: 43 86 87 88
Fax: 43 86 87 90
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: 09 80 45 51
Fax: 09 80 45 52 00
Email: info.fi@vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30 (0,18 € TTC/min)
Fax: 0 825 02 30 35 (0,18 € TTC/min)
Email: info.fr@vwr.com

Germany

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
FreeCall: 0800 702 00 07
Fax: 0180 570 22 22*
Email: info.de@vwr.com
*0,14 €/Min. aus d. dt. Festnetz

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: (52) 521-130
Fax: (52) 470-069
Email: info.hu@vwr.com

Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern
Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
Email: sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311
Fax: 800 152999/02-40090010
Email: info.it@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
Email: info.nl@vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 22 90 00 00
Fax: 815 00 940
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200
Fax: 058 32 38 205
Email: info.pl@vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Alfragide
Tel.: 21 3600 770
Fax: 21 3600 798/9
Email: info.pt@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
Email: kundservice.se@vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
Email: info.ch@vwr.com

Turkey

VWR International Laboratuvar
Teknolojileri Ltd.Şti.
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90 216 598 2900
Fax: +90 216 598 2907
Email: info.tr@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
Email: uksales@vwr.com

China

VWR International China Co., Ltd.
Shanghai Branch
Room 256, No. 3058 Pusan Road
Pudong New District
Shanghai 200123
Tel.: +86-21-5898 6888
Fax: +86-21-5855 8801
Email: info_china@vwr.com

India

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India - 560058
Tel.: +91-80-28078400
Fax: +91-80-28078410
Email: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Fax: +65 6264 3780
Email: sales.sg@vwr.com

GO TO **VWR.COM** FOR THE
LATEST NEWS, SPECIAL OFFERS
AND DETAILS OF YOUR LOCAL
VWR DISTRIBUTOR