



Product Specification Sheet

VWR DISPOSABLE WEIGHING BOATS, STANDARD DIAMOND SHAPED

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity (ml)	WxDxH (mm)	Colour	Pk
611-1993	5	55×35×6	White	500
611-1996	5	55×35×6	Black	500
611-0097	25	80×50×14	White	500
611-9184	25	80×50×14	Black	500
611-9161	30	80×60×14	White	500
611-9185	30	80×60×14	Black	500
611-9162	100	125×100×19	White	500
611-9186	100	125×100×19	Black	500

Material: **Polystyrene, Diamond shaped**

Type: **Standard**

Temperature Resistant: **Up to +70°C**

Chemical resistance: **Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols and bases. May be used with microwave moisture analyzers.**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **December 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0