



ZPACK13247010 PLASTIC BUCKET 10 L UN Y



SPECIFICATIONS

Material :	HDPE White
Filling Opening :	258,00 mm
Top Diameter :	300,00 mm
Bottom Diameter :	227,00 mm
Height :	263,00mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343