



Material	191536N
Material description	Sodium hydroxide 40 % W/V
Grade	low in nitrogen
CAS Number	1310-73-2
Molecular formula	NaOH
Molecular mass	40.00

Characteristics	Specifications
Assay (w/v)	38.6 --> 41.7 %
Appearance	Clear colourless liquid
Identification (Na)	Passes test
CO ₃ (as Na ₂ CO ₃)	Max. 0.2 %
Total N (Nitrogen)	Max. 10 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.