



Material	22610.132
Material description	Cetyltrimethylammonium bromide
Grade	AnalaR NORMAPUR
CAS Number	57-09-0
Molecular formula	CH ₃ (CH ₂) ₁₅ N(Br)(CH ₃) ₃
Molecular mass	364.45

Characteristics	Specifications
Assay	Min. 99.0 %
Heavy metals (as Pb)	Max. 10 ppm
Ignition residue (SO ₄)	Max. 0.1 %
Water	Max. 1.0 %
Fe (Iron)	Max. 10 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
VWR International bvba; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.