



Certificate of conformity

VWR DISPOSABLE WEIGHING BOATS, ANTI-STATIC VERSION DIAMOND SHAPED

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Capacity (ml)	W×D×H (mm)	Colour	Pk
611-1994	5	55×35×6	White	500
611-1995	5	55×35×6	Black	500
611-9187	25	80×50×14	White	500
611-9190	25	80×50×14	Black	500
611-9188	30	80×60×14	White	500
611-9191	30	80×60×14	Black	500
611-9189	100	125×100×19	White	500
611-9192	100	125×100×19	Black	500

Material: **Polystyrene, Diamond shaped**

Type: **Anti-static**

Temperature Resistant: **Up to +70°C**

Chemical resistance: **Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols and bases. May be used with microwave moisture analyzers.**

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2015**