



Certificate of Conformity

VWR PET WATER SAMPLING BOTTLES, SQUARE, WITH TAMPER CAP

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Capacity(ml)	Sterile	Height (mm)	Width (mm)	Cap colour / thiosulfate mg/l	PK
331-0278	500	+	192	70x70	Blue	108
331-0280	500	+	192	70x70	Red / 20mg/L	108
331-0282	500	+	192	70x70	White/ 120mg/L	108
331-0279	1000	+	222	89x89	Blue	70
331-0281	1000	+	222	89x89	Red / 20mg/L	70
331-0283	1000	+	222	89x89	White / 120mg/L	70

Material: PET for bottles and PP for caps

Properties:

Clear

Rigid

Leakproof

Do not autoclave!

Temperature Resistance: - 10 to +40 °C

Do not expose to open flame, or heaters !

Other Characteristics:

- Space saving bottles, easy to fill
- Clear PET ensures good visibility of the content level
- Ideal for water testing
- Capped and Irradiated
- Excellent shock resistance
- Dosed with Sodium Thiosulphate where required



Certificate of Conformity

Product compliance:

Sterilised by irradiation with full batch traceability

Usage temperature: Room temperature

Storage conditions: Room temperature (0°C to + 30°C) Do not expose to direct sunlight.

Autoclavable: No

Expiry: 2 years

Single use product

Product compliance:

Absence of heavy metals: Raw material is compliant with 94/62/EC directive

Food contact: Raw material is compliant to Regulation EC N°1935/2004

Standard NF en 13972: In rigid plastic containers. Nominal capacity definition, brim and total and measurement capabilities brim and total.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **December 2014**