



Product Specification



Material 350864X
Material description GIEMSA'S STAIN IMPROVED R66 SOLN.GURR

Characteristics	Specifications
Suitability for microscopy (blood smear)	Passes test
Erythrocytes	Beige-brownish
Nuclei	Violet
Eosinophilic granules	Reddish to red-brown
Neutrophilic granules	Red-violet
Lymphocyte cytoplasm	Blue
Monocyte plasma	Grey-blue
Thrombocytes	Violet
Basophilic granules	Dark violet
Blood parasite nuclei	Light red

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.