

| | |
|----------------------|-------------------|
| Material | 330884Y |
| Material description | Resazurin tablets |
| Grade | for milk testing |

Characteristics

Appearance
Shape
Average tablet weight
Diameter
Solubility in water (2 tablets/100 ml)

Specifications

Red to purple tablets
Round biconvex
240 - 260 mg
7.5 - 8.0 mm
≤ 4 min

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.