

Material	23539.291
Material description	1,4-Dioxan
Grade	AnalaR NORMAPUR for synthesis dehydrated with less than 0,0100 % of water
CAS Number	123-91-1
Molecular formula	C ₄ H ₈ O ₂
Molecular mass	88,11

Characteristics	Specifications
Assay	Min. 99,5 %
Acidity or alkalinity	Max. 0,001 meq/g
Colouration	Max. 10 APHA
Peroxides (as H ₂ O ₂)	Max. 0,0100 %
Water	Max. 0,0100 %
Acetyls and carbonyl derivatives	Max. 40 ppm
Evaporation residue	Max. 50 ppm
Fe (Iron)	Max. 1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International