

|                      |  |
|----------------------|--|
| Material             | 330902Q                                      |
| Material description | Ringer's solution tablets (quarter strength) |
| Grade                | for milk testing                             |

## Characteristics

## Specifications

|                                |                |
|--------------------------------|----------------|
| Appearance                     | White tablets  |
| Shape                          | Round biconvex |
| Clarity of solution            | Passes test    |
| Weight uniformity:             | .              |
| Max. 10 % outside 1150–1337 mg | Conforms       |
| None outside 1120–1368 mg      | Conforms       |
| Average tablet weight          | 1181 - 1306 mg |
| Diameter                       | 14.0 - 14.5 mm |
| Calcium (per tablet)           | 9.7 - 11.7 mg  |
| Chloride (per tablet)          | 700 - 750 mg   |
| Hardness                       | ≥ 4 kp         |
| pH                             | 7 - 8          |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.