



Certificate of Conformity

VWR SWIVEL CLAMPS, TALON®

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht

Cat. No.	Description	Material	Spanning Width (mm)	Length (mm)	Pk
241-0407	Dual adjustment clamp, three pronged, medium	Stainless steel	0 - 69	178	1
241-0408	Dual adjustment clamp, two pronged, medium	Stainless steel	0 - 76	163	1
241-0410	Dual adjustment clamp, three pronged, medium	Nickel-plated zinc	0 - 69	178	1

Other characteristics:

- Holds apparatus near the lab frame.
- For circular or irregular objects.
- Built in holder grips rods up to 19 mm in diameter.
- Holder is adjustable for forward or reverse facing adjustment screws.
- Can be locked in place once desired position is achieved.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2017**