



ZPACK13250250 PLASTIC BOTTLE POWDER 250 ML DIN 60



SPECIFICATIONS

Material :	HDPE White
Neck Type :	DIN 60
Height :	103,00 mm
Length :	75,00 mm
Wide :	75,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343