



Product Specification Sheet

VWR CLOSURES FOR TEST AND CULTURE TUBES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Cap size	Pk
217-7004	10 mm	3.000
217-7005	11 mm	2.500

Material: LDPE

Other Features:

Physical resistance: Flexible, high mechanical resistance

Use temperature main body: -80°C/+75°C

Autoclavage of the main item: N/A

Absence of heavy metals: Raw material is compliant with the EC/94/62 directive

Food contact: Raw material is compliant with regulation EC N°1935/2004

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **May 2014**