



| | |
|------------------------|-----------------------------------|
| Material | 30025.320 |
| Material description | Hydrochloric acid 2 mol/l (2 N) |
| Grade | AVS TITRINORM Volumetric Solution |
| CAS Number | 7647-01-0 |
| Molecular formula | HCl |
| Molecular mass | 36.46 |
| Additional information | |

| Characteristics | Specifications |
|---|-----------------------|
| Titer (20°C; real value 0.2 % accuracy) | 1.996 --> 2.004 mol/l |
| NIST traceable | Confirmed |

| Signature | Additional information |
|--|--|
| <p>We certify that this batch conforms to the specifications listed above.</p> <p>This document has been produced electronically and is valid without a signature.</p> <p>Isabelle Habay, Head of Laboratory - Briare VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare; France</p> | <p>The guaranteed value is traceable to primary Standard Reference Materials (SRM) from National Institute of Standards and Technology (NIST).</p> |

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.