



Material	83636.350
Material description	N,N-Dimethylacetamide
Grade	HiPerSolv CHROMANORM for HPLC
CAS Number	127-19-5
Molecular formula	C ₄ H ₉ NO
Molecular mass	87.12

Characteristics	Specifications
Assay (GC)	Min. 99.50 %
Evaporation residue	Max. 0.0005 %
Water	Max. 500 ppm
Transmittance (270 nm)	Min. 10.0 %
Transmittance (290 nm)	Min. 71.0 %
Transmittance (310 nm)	Min. 89.0 %
Transmittance (320 nm)	Min. 94.0 %
Transmittance (400 nm)	Min. 98.0 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International