



Certificate of Conformity

VWR FILTERING FLASK, WITH GLASS HOSE CONNECTION

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Capacity (ml)	Outer diameter 'D'±1	Inner diameter 'D2'≈	Total height 'H'±2	Thickness 'S' Min.	Height 'H2'≈	Heel radius "R" Min.	Pk
511-0260	100	70	17	107	1.5	25	6	1
511-0261	250	84	28	158	1.8	40	12	1
511-0262	500	106	33.5	188	2.0	45	15	1
511-0263	1000	135±3	37	242±3	2.0	50	20	1
511-0264	2000	165±3	41	288±4	2.5	55	35	1


Material : Borosilicate glass 3.3

Other Characteristics:

- Flasks are suitable for hot and cold filtration.
- Highly resistant against chemicals and thermal shock.
- High temperature resistance.
- With white inscription field.



Certificate of Conformity



Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2017**